AGROLABO
Canine Coronavirus Ab ELISA

Canine Coronavirus Ab ELISA

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Canine Coronavirus Ab ELISA

Kit ELISA for the detection of antibodies anti-Canine Coronavirus in canine serum or plasma samples

TEST PRINCIPLE

This test is performed as an indirect enzyme linked immunosorbent assay (indirect ELISA). The antigen of Canine Coronavirus is fixed on the wells of the plate. The samples to be analyzed are distributed in each well. If the sample contains specific antibodies anti-Canine Coronavirus, these antibodies will bind to the antigen adsorbed, forming the Ag-Ab complex. After washing to eliminate all non-fixed material, the presence of anti-Canine Coronavirus antibodies is detected by specific anti-canine IgG monoclonal antibody, labelled with horseradish peroxidase (conjugate). All non-fixed material is removed by a second washing. A specific enzyme substrate is then added to the wells which causes a colorimetric reaction when the conjugate is present.

TEST COMPONENTS

- 1 x 96 wells plate (in strip of 8 wells);
- 1 dropper of positive control;
- 1 dropper of negative control;
- 1 dropper of conjugate - 100X;
- 1 bottle of washing solution - 10X;
- 1 bottle of diluent - ready for use;
- 1 dropper of substrate (TMB) - ready for use;
- 1 dropper of stop solution - ready for use.

MATERIALS REQUIRED BUT NOT SUPPLIED

Test tubes for sample dilutions; precision pipettes and tips; spectrophotometer at a wavelength of 450 nm.

PRECAUTIONS

1. Read the instructions carefully.
2. Bring all reagents to room temperature (+20-25°C) prior to use.
3. Do not substitute or mix instructions or reagents from different kits or lots.
4. The substrate is very sensitive to the light and to contaminations. Do not use the tip directly in the vial. The stop solution is a diluted strong acid. Handle the two reagents with care.
5. Do not use components after expiry date.
6. Do not smoke, eat or drink where specimens or kit reagents are being handled.
7. In each utilization of the kit the positive and negative controls must be tested.
8. Use a new tip for each sample.
9. Observe each incubation time.
10. Prepare each sample according to instructions.
11. Avoid any contamination of the reagents.
12. All kit components must be stored at +2-8°C.

WASHING STEPS

The washing steps may be carried out using an automated plate washer or a multichannel pipetting device suitable for dispensing 300 µl in each well. After the incubation, the washing steps should be carried out according to the following instructions:

- throw out the content of the plate by vigorously turning the plate over to avoid the possible mixture of the content from one well to another,
- dispense a volume of 300 µl of wash buffer into each well,
- shake the plate gently, avoiding contamination between wells,
- turn over the plate and tap the wells vigorously to remove the wash buffer,
- repeat the process as many times as is indicated by the kit's instructions,
- prior to emptying out the plate's contents during the final washing step, ensure that the next reagent to be added to the plate is ready to use,
- the wells should not be left dry for any longer than strictly necessary in order to dispense the reagent of the following step,
- after the last washing step, blot the empty wells face down on clean paper towels or absorbent paper to remove residual wash buffer.

PREPARATION OF SAMPLES

Canine serum or plasma samples.

For the titration, we recommend to make at least four two-fold dilutions for each sample (1:50; 1:100; 1:200; 1:400).

PREPARATION OF REAGENTS

Wash buffer solution:

Dilute one part of the wash buffer supplied into 9 parts of distilled or deionized water (ie. 100 ml of concentrated wash buffer with 900 ml of water). Once diluted, the wash buffer solution may be stored at +2-8°C.

Conjugate:

To be prepared immediately prior to use. Dilute the required amount of conjugate **1:100** with the diluent. Mix the solution prior to use. Prepare only the required amount and do not reuse surplus conjugate.

Positive and Negative Controls:

Dilute as samples (1:50; 1:100; 1:200; 1:400).

TEST PROCEDURE

1. Bring all reagents to room temperature before use.
2. Add 100 µl of each positive control dilutions to the first 4 wells and 100 µl of each negative control dilutions to the following 4 wells. Add 100 µl of each sample dilutions to the remainder wells of the plate.

Seal the plate and incubate at room temperature for 10 minutes.

3. Wash cycle 3 times or more according to the procedure previously described.
4. Dispense 100 µl of previously prepared conjugate into each well. Delicately shake the plate to obtain a correct homogenization of the reagents.

Seal the plate and incubate at room temperature for 10 minutes.

5. Wash cycle 5 times or more according to the procedure previously described.
6. Add 100 µl of substrate to each well

Seal the plate and incubate at room temperature for 5 minutes.

7. Add 100 µl of stop solution to each well.
8. Measure the absorbance value (OD) of each well with a spectrophotometer at a wavelength of 450 nm.

READING AND INTERPRETATION OF THE RESULTS

The reading of the absorbance values (OD) must be carried out at a wavelength of 450 nm. If you run the samples or controls in duplicate wells, the OD value for the sample and controls will be the mean of both wells.

Test Validation

The test is to be considered valid when:

Titre of Positive control: OD \geq 1:200

OD Negative control: OD < 0.4

CUT-OFF value determination

The cut-off value is OD = 0.4

Interpretation of the results

Positive samples: samples with an OD higher than 0.4

Negative samples: samples with an OD lower than 0.4

TITRATION

The antibody titre is represented by the highest dilution of the sample that has an OD value higher than 0.4

AGROLABO

Manufactured by;

Agrolabo SpA

Head Office, Laboratories

and Production Centre

Diagnostic Division

Via Masero 59

10100 Scarmagno (TO)

Italy

Tel +39 0125 731111

Fax +39 0125 731190

E-mail agrolabo@agrolabo.it

www.agrolabo.it

shop.agrolabo.it