Professional Veterinary Diagnostic Test  
**BioPronix** product line  
 For Veterinary use only  
 For in vitro use only  
 Manufactured by Agrolabo S.p.A.

## INTERPRETATION OF THE RESULTS

### Negative Result

The test is negative if a red line appears in window C (control line).

### Positive Result

The test is positive if two red lines appear. One line will appear in window T (test line) and one in window C (control line).

### Invalid Result

The test is invalid if no coloured band appears in window C, even if a coloured band appears in window T.


## QUICK GUIDE

1. Take a sample and place it in the tube with the buffer solution. Shake carefully then place the tube in a vertical position.
2. Once the residues have settled, collect the supernatant with one of the pipettes provided in the kit and add 3-4 drops into the S window.
3. Wait 10 minutes to read the results, not later than 15 minutes.

4. Interpretation of the results.

## WARNING


Open the test package immediately before use. Place the device horizontally for testing. Keep the pipette in a vertical position when adding the sample to the device.

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
Invalid

## DISTEMPER IC

### Rapid test kit for the detection of Canine Distemper Virus

## USER MANUAL

Rev. R1 – 16/02/21

### TEST PRINCIPLE

Distemper IC test employs Immunochromatographic sandwich technique using two monoclonal antibodies against the antigens of Canine Distemper Virus (CDV). A monoclonal antibody is conjugated with colloidal gold and a further monoclonal antibody is immobilized on the membrane. If the sample contains CDV antigens (positive sample), these antigens will bind to the colloidal gold conjugated antibody forming the Ab-Ag complex. The complex migrates along the membrane and will be captured by the second monoclonal antibody immobilized on the membrane at T window level, forming a red line (test line). The sample continues to migrate along the membrane up to window C, where a reagent able to bind the antibody conjugated with colloidal gold is immobilized, resulting in the formation of a red line (control line).

### KIT COMPONENTS

- Disposable pipettes
- Tubes with diluent
- Sterile swabs
- Labels for tubes identification
- DISTEMPER IC Devices

### STORAGE AND SHELF LIFE

Store at room temperature (+18-25°C).

When properly stored, the product remains stable for its full shelf life.

Keep the kit away from excessive heat/cold or direct sunlight.

### PRECAUTIONS

IMPORTANT: open individual device packaging immediately prior to use.

For veterinary use only.

Do not use kit components after expiry date.

The test device should be placed horizontally on a flat surface while testing.

Use a new disposable pipette for each sample.

### SAMPLE

Canine conjunctiva or snivel: take the sample with a sterile swab and put into a tube with the diluent.

The samples can be stored at +2-8 °C for no longer than 7 days.

All samples need to be brought at room temperature before use.

### TEST PROCEDURE

1. Remove the test device from its individual wrapping, placing it on a flat surface.

2. Take a sample dipping the sterile swab. Insert the sample into the appropriate tube containing diluent. Carefully shake the tube to mix the substances and obtain the correct sample dilution. Place the tube in a vertical position, in order to allow residues to precipitate.

3. Once the residues have precipitated, use a new pipette to extract the supernatant from the tube. Holding the pipette approximately 1 cm from the test device, add 3-4 drops of supernatant to window S of the device.

4. Read the results after 10 minutes and no later than 15 minutes.