AGROLABO
DISTEMPER IgG Ab ELISA

DISTEMPER IgG Ab ELISA

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DISTEMPER IgG Ab ELISA

Kit ELISA for the detection of IgG antibodies anti-Canine Distemper virus in canine serum or plasma samples

TEST PRINCIPLE

This test is performed as an indirect enzyme linked immunosorbent assay (indirect ELISA). The antigen of Canine Distemper Virus (CDV) is fixed on the wells of the plate. The samples to be analyzed are distributed in each well. If the sample contains specific antibodies anti-CDV, these antibodies will bind to the antigen adsorbed, forming the Ag-Ab complex. After washing to eliminate all non-fixed material, the presence of anti-CDV antibodies is detected by a specific anti-canine IgG monoclonal antibody, labelled with horseradish peroxidase (conjugate). All non-fixed material is removed by a second washing. A specific enzyme substrate is then added to the wells which causes a colorimetric reaction when the conjugate is present.

TEST COMPONENTS

- 1 x 96 wells plate (in strip of 8 wells);
- 1 dropper of positive control - ready for use;
- 1 dropper of negative control - ready for use;
- 1 dropper of conjugate - ready for use;
- 1 bottle of washing solution - 10X;
- 1 bottle of diluent - ready for use;
- 1 dropper of substrate (TMB) - ready for use;
- 1 dropper of stop solution - ready for use.

MATERIALS REQUIRED BUT NOT SUPPLIED

Test tubes for sample dilutions; precision pipettes and tips; spectrophotometer at a wavelength of 450 nm.

PRECAUTIONS

1. Read the instructions carefully.
2. Bring all reagents to room temperature (+20-25°C) prior to use.
3. Do not substitute or mix instructions or reagents from different kits or lots.
4. The substrate is very sensitive to the light and to contaminations. Do not use the tip directly in the vial. The stop solution is a diluted strong acid. Handle the two reagents with care.
5. Do not use components after expiry date.
6. Do not smoke, eat or drink where specimens or kit reagents are being handled.
7. In each utilization of the kit the positive and negative controls must be tested.
8. Use a new tip for each sample.
9. Observe each incubation time.
10. Prepare each sample according to instructions.
11. Avoid any contamination of the reagents.
12. All kit components must be stored at +2-8°C.

WASHING STEPS

The washing steps may be carried out using an automated plate washer or a multichannel pipetting device suitable for dispensing 300 µl in each well. After the incubation, the washing steps should be carried out according to the following instructions:

- throw out the content of the plate by vigorously turning the plate over to avoid the possible mixture of the content from one well to another,
- dispense a volume of 300 µl of wash buffer into each well,
- shake the plate gently, avoiding contamination between wells,
- turn over the plate and tap the wells vigorously to remove the wash buffer,
- repeat the process as many times as is indicated by the kit's instructions,
- prior to emptying out the plate's contents during the final washing step, ensure that the next reagent to be added to the plate is ready to use,
- the wells should not be left dry for any longer than strictly necessary in order to dispense the reagent of the following step,
- after the last washing step, blot the empty wells face down on clean paper towels or absorbent paper to remove residual wash buffer.

PREPARATION OF SAMPLES

Canine serum or plasma samples.

Prepare a **1:100** dilution of each sample to be tested using diluent.

PREPARATION OF REAGENTS

Wash buffer solution:

Dilute one part of the wash buffer supplied into 9 parts of distilled or deionized water (ie. 100 ml of concentrated wash buffer with 900 ml of water). Once diluted, the wash buffer solution may be stored at +2-8°C.

TEST PROCEDURE

1. Bring all reagents to room temperature before use.
2. Add 100 µl of positive and negative controls in two different wells.
Add 100 µl of diluted samples on the remainder wells of the plate.
Seal the plate and incubate at room temperature for 10 minutes.
3. Wash cycle 4 times or more according to the procedure previously described.
4. Add 100 µl of conjugate to each well. Delicately shake the plate to obtain a correct homogenization of the reagents.
Seal the plate and incubate at room temperature for 10 minutes.
5. Repeat the wash cycle 4 times or more according to the procedure previously described.
6. Add 100 µl of substrate to each well.
Seal the plate and incubate at room temperature for 5 minutes.
7. Add 100 µl of stop solution to each well.
8. Measure the absorbance value of each well with a spectrophotometer at a wavelength of 450 nm.

READING AND INTERPRETATION OF THE RESULTS

The reading of the absorbance values (OD) must be carried out at a wavelength of 450 nm. If you run the samples or controls in duplicate wells, the OD value for the samples and controls will be the mean of both wells.

Test Validation:

The test is to be considered valid when:

Positive Control: $OD > 1$

Negative Control: $OD < 0.2$

CUT-OFF value determination:

Cut-off = OD positive control x 0.2

INTERPRETATION OF THE RESULTS

Positive samples: samples with an OD higher than the cut-off value.

Negative samples: samples with an OD lower than the cut-off value.

TITRATION:

Low titers (corresponding to IFA values of 1:20-1:40): samples with OD values ranging from 0.2 x OD positive control to 0.4 x OD positive control.

Medium titers (corresponding to IFA values of 1:80-1:160): samples with OD values ranging from 0.4 x OD positive control to 0.8 x OD positive control.

High titers (corresponding to IFA values $\geq 1:320$): samples with OD values higher than 0.8xOD positive control.

AGROLABO

Manufactured by;

Agrolabo SpA

Head Office, Laboratories

and Production Centre

Diagnostic Division

Via Masero 59

10100 Scarmagno (TO)

Italy

Tel +39 0125 731111

Fax +39 0125 731190

E-mail agrolabo@agrolabo.it

www.agrolabo.it

shop.agrolabo.it