INTERPRETATION OF THE RESULTS

Negative Results:

The test is negative if a red line appears in windows C (control line).

Positive Results Toxo IgM:

The test is positive for Toxo IgM if two red lines appear in column 1. One line will appear in window T (test line) and one in window C (control line).

Positive Results Toxo IgG:

The test is positive for Toxo IgG if two red lines appear in column 2. One line will appear in window T (test line) and one in window C (control line).

Invalid Results:

The test is invalid if no red band appears in window C, even if a red band appears in window T.


QUICK GUIDE

1. Remove the device from its wrapping and place it on a flat surface. Use the micropipette to add 10 µl of serum, plasma or whole blood with anticoagulant into both windows S of the device.

2. Holding the diluent bottle vertically add 2 drops of diluent into both windows S of the device.

3. Wait 5-10 minutes to read the results, not later than 15 minutes.

4. Interpretation of the results.

WARNING

Open the test package immediately before use. Place the device horizontally for testing. Keep the pipette in a vertical position when adding the sample to the device.

Professional Veterinary Diagnostic Test
BioPronix product line
For Veterinary use only
For in vitro use only
Manufactured by Agrolabo S.p.A.

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TEST PRINCIPLE

Feline Toxo IgM-IgG IC is an immunochromatographic test for the detection of anti-Toxoplasma gondii IgM and IgG antibodies in feline whole blood, serum or plasma.

If the sample contains anti-Toxoplasma IgM and IgG antibodies (positive sample), it migrates to the T window, where these antibodies will bind to the Toxoplasma Ag immobilized on the membrane, forming the Ab-Ag complex. The addition of the diluent will release the antibody (Ab) anti-cat IgG or anti-cat IgM conjugated with colloidal gold that will specifically bind the Ab-Ag complex present on the membrane, forming a red line in window T.

The sample continues to migrate along the membrane up to window C, where a reagent able to bind the antibody conjugated with colloidal gold is immobilized resulting in the formation of a red line (control line).

KIT COMPONENTS

- 10 µl Micropipette
- Tips
- Vial with diluent
- FELINE TOXO IgM-IgG IC devices

STORAGE AND SHELF LIFE

Store at room temperature (+18-25°C). When properly stored, the product remains stable for its full shelf life. Keep the kit away from excessive heat/cold or direct sunlight.

PRECAUTIONS

IMPORTANT: open individual device packaging immediately prior to use.

For veterinary use only.

Do not use kit components after expiry date.

The test device should be placed horizontally on a flat surface while testing.

Use a new tip for each sample.

SAMPLES

All samples need to be brought to room temperature before use. Previously frozen or stored samples must be centrifuged before use.

Feline whole blood samples with anticoagulant (EDTA, heparin, etc) may be used. Whole blood samples should be fresh or stored in a refrigerator for no longer than 5 days.

Serum or plasma that is fresh, or that has been refrigerated (+2-8°C) for no longer than 5 days, can be used. For longer storage, samples should be frozen.

Haemolyzed samples may affect test results.

TEST PROCEDURE

1. Remove the device from its individual wrapping and place it on a flat surface.

2. Holding the micropipette add 10 µl of whole blood with anticoagulant, serum or plasma into both windows S of the device.

3. Holding the diluent bottle in a vertical position add 2 drops into both windows S of the device.

4. Read the result after 5-10 minutes and no later than 15 minutes.