# FLUO ANA HEp-2

Kit IFA for the detection of anti-nuclear antibodies  
in canine serum or plasma

## INSTRUCTIONS FOR USE

### INTENDED USE

The FLUO ANA HEp-2 kit is intended for the semiquantitative detection of anti-nuclear antibodies in canine serum or plasma by indirect immunofluorescence assay.

### TEST PRINCIPLE

Teflon-masked substrate slides are coated with a specific HEp-2 cell culture. Canine sera are diluted in buffered saline and incubated in the individual slide wells to allow reaction of patient anti-nuclear antibody (ANA) with the cell culture coated on the wells. Slides are then washed to remove unreacted serum proteins, a fluorescence-labelled FITC anti-canine IgG (conjugate) is added. This conjugate is allowed time to react with antigen-antibody complexes. The slides are washed again to remove unreacted conjugate. The resulting reactions can be visualized using standard fluorescence microscopy, where a positive reaction is seen as sharply defined apple-green fluorescent pattern described after in the paragraph of test results.

A negative reaction is seen with fluorescence unlike that seen in the positive control well (red-counterstained cells). Positive reactions may then be retested at a higher dilution to determine the highest reactive (endpoint titer).

### TEST COMPONENTS

- 10 x 10 wells substrate slides;
- 1 vial of FITC anti-dog IgG conjugate-ready for use;
- 1 vial of positive control-ready for use;
- 1 vial of negative control- ready for use;
- 1 vial of mounting medium-ready for use;

### MATERIALS REQUIRED BUT NOT SUPPLIED

PBS buffer (phosphate buffered saline 1X pH 7.2-7.4); test tubes for serum dilutions; precision pipettes and tips; 24x50 mm coverslips; fluorescence microscope with filter system for FITC (fluorescein isothiocyanate, excitation wavelength 465-495, barrier filter 515-555) and 400X magnification; +37°C incubator; humid chamber for slides.

### STORAGE

Kit components should be stored at +2-8°C. Bring them to room temperature (+20-25 °C) before performing the test.

### SAMPLE

Canine serum or plasma. Allow sample to clot and separate serum by

centrifugation. Transfer serum aseptically to a tightly closing sterile container. Store at +2-8 °C. If testing is to be delayed longer than 5 days, freezing the sample at -20 °C or colder is recommended. Acute specimens should be drawn at the onset of illness; convalescent specimens should be obtained at two and four weeks intervals to check for titer changes.


### TEST PROCEDURE

- Prepare **1:50** screening dilution for all untested sera. For sera found positive on a previous assay run, prepare serial two-fold dilution in PBS, starting from the dilution 1:50.
- The positive and negative controls are ready to use. Do not dilute them. For each serum dilution to be tested add 20 µl to one slide well and record the location for later reference. For each assay run include the negative control and the positive control (20 µl).
- Place slides in a humid chamber and incubate for 30 min at +37°C.
- Remove humid chamber from incubator.
- Washing step: Tap remaining serum dilutions gently from the slides and shake the slides gently for 5 minutes in PBS. Repeat this step for another 5 minutes with fresh PBS. Briefly rinse the slides with distilled water. Tap remaining water gently from the slides and, if necessary, dry the Teflon mask between the wells with absorbent paper or cotton bud sticks. However, do not allow the antigen wells to dry out. If using a washing bottle, do not focus the stream directly onto the antigen wells.
- Add to each slide well 1 drop (20 µl) of conjugate, then return slides to the humid chamber for another 30 minutes incubation at +37°C. Incubation should be in the dark to protect the photosensitive conjugate.
- Repeat washing steps as previously described.
- Add 2 drops of mounting fluid to each slide and place the coverslip on, carefully removing air bubbles caught under the coverslip.
- Evaluate the slides at 400X magnification, comparing each well to the visual intensity and appearance of the Leishmania seen in the positive and negative control wells. Slides may be stored at +2-8°C in the dark for up to 7 days.

### RESULTS

For the evaluation, a fluorescence microscope with a filter system for FITC and 400X magnification is required. The fluorescence pattern (form, density etc.) of the Negative and Positive Control is considered as reference pattern. Patterns of reactivity different than those seen in the controls must be considered non-specific, which means negative.

**Positive test result at  $\geq$  1:50 screening dilution:** a positive result could show different patterns of nuclear and/or cytoplasmic apple green fluorescence, but any patterns observed at 1:50 dilution or higher

**Negative test result at 1:50 screening dilution:** Report as negative for antinuclear antibody. No yellow-green fluorescence is visible.

**Homogeneous pattern:** seen as homogeneous or diffuse green fluorescence of the nucleus due to the autoantibodies to native DNA histories and/or DNP (deoxyribonucleoprotein). The chromosome of the mitotic cell (dividing cells) are important indicators of a homogeneous pattern because they will stain as irregularly shaped masses with more intensely stained outeredges.

**Speckled pattern (most commonly):** cells not in division (non-mitotic) will show "true speckled" green fluorescence patterns seen with centromere antibodies. The chromosomal area of mitotic cells will be not stained, the cytoplasm of them shows a few speckles.

**Nucleolar pattern:** seen as homogeneous, clumped or speckled staining of the nucleolus, often associated with a dull, homogeneous fluorescence in the rest of the nucleolus. In the mitotic cells the chromosome will be negative.


**Peripheral pattern:** seen as predominant staining of the nucleoli periphery. A further dilution of the positive samples is recommended to determine the highest dilution that is still positive (endpoint titer). An acute infection (2 to 4-fold titer increase) can only be determined by the titer determination of a coupled serum test (2 samples in an interval of 2-3 weeks). The interpretation of test results should always be based on anamnestic and especially clinical data and additional laboratory parameters.

### QUALITY CONTROL

The negative control serum and the positive control serum should be assayed with each daily run. The negative control well is an example of a non-reactive serum, with either uniform red counterstain or slight, but uniform greenish staining. The fluorescence intensity of the positive control is an example of medium intensity fluorescence. If one of the controls does not react as specified, the assay run should be considered invalid, reagent components and procedural steps should be re-checked, and the assay repeated.

### PRECAUTIONS

For veterinary use only. Follow the instructions carefully. Since no testing can assure the absence of infectious agents, these reagents, as well as all serum specimens and equipment coming in contact with these specimens, should be handled with good laboratory practices to avoid skin contact and ingestion. The substrate slides are prepared with chemically inactivated antigens. However, the slides should be considered potentially infectious and handled accordingly. Conjugate is photosensitive, for protection store in the dark and return to storage after use. Conjugate contains Evans blue dye, avoid ingestion and skin contact. Do not use after the expiration date.

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