Professional Veterinary Diagnostic Test

For Veterinary use only
For in vitro use only

Manufactured by ACCUPLEX DIAGNOSTICS Ltd.
Distributed exclusively in Italy by AGROLABO SpA.


INTERPRETATION OF THE RESULTS

Note: the intensity of the colour on each line may vary depending on the level of SAA in the blood sample.

Horses with normal levels of protein SAA will give 4 lines in the test window. As the severity of inflammation increases so will the level of protein SAA. The test lines will get progressively weaker in colour so that only 3 lines become visible indicating that the horse has mild inflammation, while 2 lines indicating that the horse has moderate inflammation. With significant inflammation, only 1 line will be visible. If no line appears the test is invalid, repeat the test using a new device.

QUICK GUIDE

1. Before starting the test remove test cassette from foil pouch
2. Touch sample applicator to blood in tube or at puncture site. DO NOT immerse sample applicator in blood. The correct volume of blood will be drawn up into the cup
3. Touch sample applicator to the sample window with a little pressure to dispel blood
4. Holding the diluent bottle in a vertical position, add 3 drops of diluent into the sample window
5. Read the result after 10 minutes, not later than 15 minutes
6. Interpretation of the results

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
EQUICHEK™-SAA - Rapid test for the detection of inflammation in equine blood

USER MANUAL

Rev. R2-27/04/22

TEST PRINCIPLE

Serum amyloid A (SAA) is a major acute phase protein of inflammatory processes in horses. It is present at very low levels when the horse is in normal health conditions, but, within a few hours of the onset of any pathogenic event, it can increase up to 10-100 times above normal. The increase in SAA during an inflammatory response is unequivocal, high and rapid. Monitoring the level of the SAA protein has a very significant advantage, since, once the treatment picture has been successfully concluded, its level rapidly decreases returning to normal values. The lowering of the SAA level therefore represents a clear sign of recovery. EQUICHEK is a competitive test for detecting SAA in whole blood. The test provides a semi-quantitative visual result that indicates whether the horse has no inflammation or if it has weak, moderate or clinically significant inflammation. The test is based on the competition between the SAA present in the sample and the SAA adhered to the test membrane that will bind to the anti-SAA antibody labelled with colloidal gold.

KIT COMPONENTS

- Disposable capillary pipettes
- Vial with diluent
- EQUICHEK devices

STORE AND SHELF LIFE

Store at room temperature (+18-25°C). When properly stored, the product remains stable for its full shelf life. Keep the kit away from excessive heat/cold or direct sunlight.

PRECAUTIONS

1. As with all diagnostic tests, a definitive clinical diagnosis should not be based on a single result, but should only be made after all clinical

- and laboratory findings have been considered.
2. This product is intended only horse blood.
 3. Use only the diluent provided with the kit.
 4. All specimens should be considered potentially hazardous and handled as an infectious agent.
 5. Do not use this product after the expiration date indicated on the outer box.
 6. Do not use the test strip if foil pouch is damaged or already opened.

SAMPLES

The test is designed for use with freshly taken whole blood, but serum or plasma can also be used. For animal side application, a needle (no syringe attached) should be inserted into a vein. Alternatively, blood can be collected into a tube (EDTA or Heparin). Blood samples with anticoagulant can be stored at +2-8°C for 5 days before the test is performed. Do not use haemolytic samples. Serum or plasma can be used, in this case the sample quantity must be taken with a precision pipette. The volume of serum or plasma to be tested is 3 µl.

TEST PROCEDURE

1. Remove the device from its individual wrapping and place it on a flat surface.
2. Touch sample applicator to blood in tube or at puncture site. DO NOT immerse sample applicator in blood. The correct volume of blood will be drawn up into the cup.
3. Touch sample applicator to the sample window with a little pressure to dispel blood.
4. Holding the diluent bottle in a vertical position, add 3 drops of diluent into the sample window.
5. Read the result after 10 minutes and no later than 15 minutes. Observe the number of lines that appear in the test window. Results will be visible in 2-3 minutes but it is recommended to read them after 10 minutes.