# FLUO HERPESVIRUS 1 and 4

Kit IFA for the detection of anti-Equine Herpesvirus 1 and 4 IgG antibodies in equine serum or plasma

## INSTRUCTIONS FOR USE

### INTENDED USE

The FLUO HERPESVIRUS 1 and 4 kit is intended for the semiquantitative detection of IgG antibodies to Equine Herpesvirus 1 (EHV-1) and 4 (EHV-4) in equine serum or plasma by indirect immunofluorescence assay

### TEST PRINCIPLE

Substrate slides consist of Teflon-masked wells containing cells infected by EHV (in the upper line by EHV-1 and in the lower line by EHV-4). Equine sera are diluted in buffered saline and incubated in the individual slide wells to allow reaction of patient antibody with the EHV antigens on the wells.

Slides are then washed to remove unreacted serum proteins, a fluorescence-labelled FITC anti-equine IgG (conjugate) is added. This conjugate is allowed time to react with antigen-antibody complexes. The slides are washed again to remove unreacted conjugate.

The resulting reactions can be visualized using standard fluorescence microscopy, where a positive reaction is seen as sharply defined apple-green fluorescence binds to the superficial proteins in each infected cells. A negative reaction is seen as a red-counterstained cells or as a fluorescence unlike that seen in the positive control well.

Positive reactions may then be retested at a higher dilution to determine the highest reactive (endpoint titer).

### TEST COMPONENTS

- 10 x 10 wells substrate slides;
- 1 vial of FITC anti-horse IgG conjugate-ready for use;
- 1 vial of positive control-ready for use;
- 1 vial of negative control- ready for use;
- 1 vial of mounting medium-ready for use;

### MATERIALS REQUIRED BUT NOT SUPPLIED

PBS buffer (phosphate buffered saline 1X pH 7,2-7,4); test tubes for serum dilutions; precision pipettes and tips; 24x50 mm coverslips; fluorescence microscope with filter system for FITC (fluorescein isothiocyanate, excitation wavelength 465-495, barrier filter 515-555) and 400X magnification; +37°C incubator; humid chamber for slides.

### STORAGE

Kit components should be stored at +2-8°C. Bring them to room temperature (+20-25 °C) before performing the test.


### SAMPLE

Equine serum or plasma.

Allow sample to clot and separate serum by centrifugation. Transfer serum aseptically to a tightly closing sterile container. Store at +2-8 °C. If testing is to be delayed longer than 5 days, freezing the sample at -20 °C or colder is recommended. Acute specimens should be drawn at the onset of illness; convalescent specimens should be obtained at two and four weeks intervals to check for titer changes.

### TEST PROCEDURE

- Prepare different screening dilutions for all untested sera (**1:50** dilution for EHV-1 and **1:100** dilution for EHV-4). For sera found positive on a previous assay run, prepare serial two-fold dilution in PBS, starting from 1:50 or 1:100.
- The positive and negative controls are ready to use. Do not dilute them. For each serum dilution to be tested add 20 µl to one slide well EHV-1 and to one slide well EHV-4 and record the location for later reference. For each assay run include the negative control and the positive control (20 µl) in both lines of wells (EHV1- and EHV-4).
- Place slides in a humid chamber and incubate for 30 minutes at +37 °C.
- Remove humid chamber from incubator.
- Washing step: Tap remaining serum dilutions gently from the slides and shake the slides gently for 5 minutes in PBS. Repeat this step for another 5 minutes with fresh PBS. Briefly rinse the slides with distilled water. Tap remaining water gently from the slides and, if necessary, dry the Teflon mask between the wells with absorbent paper or cotton bud sticks. However, do not allow the antigen wells to dry out. If using a washing bottle, do not focus the stream directly onto the antigen wells.
- Add to each slide well 1 drop (20 µl) of conjugate, then return slides to the humid chamber for another 30 minutes incubation at +37 °C. Incubation should be in the dark to protect the photosensitive conjugate.
- Repeat washing steps as previously described.
- Add 2 drops of mounting fluid to each slide and place the coverslip on, carefully removing air bubbles caught under the coverslip.
- Evaluate the slides at 400X magnification, comparing each well to the visual intensity and appearance of the fluorescence seen in the positive and negative control wells. Slides may be stored at +2-8°C in the dark for up to 7 days.

### RESULTS

For the evaluation, a fluorescence microscope with a filter system for FITC and 400X magnification is required. The fluorescence pattern (form, density etc.) of the Negative and Positive Control is considered as reference pattern. Patterns of reactivity different than those seen in the controls must be considered non-specific, which means negative.

**Negative test result at 1:50 screening dilution x EHV-1:** report as negative for EHV-1 IgG antibodies.

**Negative test result at 1:100 screening dilution x EHV-4:** report as negative for EHV-4 IgG antibodies.

No clear visible fluorescence in the cytoplasm of the infected cells.

**Positive test result EHV-1 at ≥ 1:50 screening dilution**

**Positive test result EHV-4 at ≥ 1:100 screening dilution**

in the cytoplasm of the infected cells it is seen a clear, yellow-green fluorescence.

A further dilution of the positive samples is recommended to determine the highest dilution that is still positive (endpoint titer).

An acute infection (2 to 4-fold titer increase) can only be determined by the titer determination of a coupled serum test (2 samples in an interval of 2-3 weeks).

The prevalence of EHV-1 and EHV-4 IgG antibodies varies depending upon the geographic region and population being tested. Serological test results should be always interpreted in conjunction with case history and clinical symptoms.


### QUALITY CONTROL

The negative control serum and the positive control serum should be assayed with each daily run. The negative control well is an example of a non-reactive serum, with either uniform red counterstain or slight, but uniform greenish staining. The fluorescence intensity of the positive control is an example of medium intensity fluorescence. If one of the controls does not react as specified, the assay run should be considered invalid, reagent components and procedural steps should be re-checked, and the assay repeated.

### PRECAUTIONS

For veterinary use only. Follow the instructions carefully. Since no testing can assure the absence of infectious agents, these reagents, as well as all serum specimens and equipment coming in contact with these specimens, should be handled with good laboratory practices to avoid skin contact and ingestion. The substrate slides are prepared with chemically inactivated antigens. However, the slides should be considered potentially infectious and handled accordingly. Conjugate is photosensitive, for protection store in the dark and return to storage after use. Conjugate contains Evans blue dye, avoid ingestion and skin contact.

Do not use after the expiration date.

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