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## I. TECHNICAL BASIS

This kit is based on an indirect enzymatic immunoassay (Indirect ELISA). We make a brief description of the technique below:

We fix the antigen on a solid support (polystyrene plate). When a sample serum contains specific antibodies against the virus, they will bind to the antigen adsorbed on plate. After washing to eliminate all non-fixed material from the sera sample, we can detect the presence of swine immunoglobulins using a specific peroxidase conjugate. After addition of the substrate a colorimetric reaction will appear which could be measured by a spectrophotometer.

In this way the presence of colour means the presence of antibodies against the virus in the swine sera, and the absence of colour the absence specific antibodies.

The aim of this kit is provide to the users with a reliable and automatizable diagnostic technique for this disease, that till this time is been made by the Hemagglutination Inhibition test. In our kit is important to remark the use of antigen obtained by the expression of the viral proteins in baculovirus growth in insect cell cultures. This method warranties the total absence of infectivity in the kit.

## II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the use instructions carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the reagents of the Kit.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handle.
7. Do not pipette by mouth.
8. Use a new tip for each serum sample.
9. For each utilisation of the Kit use a fresh substrate preparation.
10. For each utilisation of the Kit, positive and negative control sera must be tested in a systematic way.
11. Substrate and substrate buffer must be handle whit care.

## III. STORAGE OF COMPONENTS

All reagents and components must be stored at between +2°C and +8°C. WARNING: Stop Solution may precipitate at this temperature, if it happens, warm up at +37°C before use.

## IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichanel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

## V. PREPARATION OF SAMPLE

- **For SCREENING POSITIVE / NEGATIVE:** Dilute the sample 1/100 in the diluent supplied.
- **For APROXIMATION TO TITRE:** Use a single dilution of 1/200 (See section IX for single-well titration approach).

## INGEZIM PPV 11.PPV.K1

- **For real TITRATION** of pig serum, we recommend to assay two-fold dilution starting from 1/100 (5µl of sample and 495 of diluent), ranging to 1/12800.
- To ease the process, this step could be done directly on the plate as follows:
- Add 200 µl of 1/100 dilution to the first well of each column.
  - Add 100 µl of diluent to the seven remaining wells of each row.
  - Transfer 100 µl from the first well to the second, mix well and proceed in the same way till the last well (discard 100µl from the last well).

## VI. PREPARATION OF REAGENTS

- **Washing solution:**  
Dilute one part of the concentrate washing solution provided in the kit with 24 parts of distilled or deionized water (i.e. 100 ml of concentrate with 2,4 L of water). When ready this solution remains stable between +2°C and +8°C.
  - **Preparation of Diluent:**  
Supplied diluent is 5x concentrated. Before using it, a 1/5 dilution must be done with distilled or deionized water (1 vol. of concentrated supplied Diluent with 4 volume of distilled water).
  - **Control sera:**  
Controls are supplied ready to use. Do not dilute
  - **Preparation of the conjugate: to make immediately before use.**  
Dilute the needed quantity of conjugate provided in the kit 1/100 into diluent:  
→ The necessary and sufficient quantity for a complete plate is 110 µl of conjugate with 10 ml of diluent.  
→ The necessary and sufficient quantity for an 8 well strip is 10µl of conjugate with 1 ml of diluent.
  - **Preparation of the substrate:**  
Ready to use.
- Shake very well the solution before the use. Prepare only the quantity needed for each time because the remainder volume has to be rejected.

## VII. TEST PROCEDURE

1. All reagents (except conjugate) must be allowed to come to room temperature before use.
2. Add 100 µl of positive control serum, to two wells of the first row of the plate and 100 µl of the negative control serum to other two wells. Add 100 µl of each dilution of sera samples to test, prepared following the previous instructions, on the remainder wells of the plate (use one row per serum). Seal the plate **and incubate for 1 h at 37°C**.
3. Wash 3 times following the procedure previously described.
4. Add 100 µl of conjugate, prepared following the previous instructions, to each well. Seal the plate and **incubate for 1 hour at 25°C** (room temperature).
5. Wash 4 times following the described procedure.
6. Add 100 µl of substrate to each well. Keep the plate **for 10 min at room temperature**.
7. Add 100 µl of stop solution to each well.
8. Read the absorbance of each well with a spectrophotometer at 405 nm within 5 min after the addition of stop solution.

## VIII. READING AND RESULT INTERPRETATION

### ◆ Validation of the plate:

Test is only valid when:

- OD Control (+) is higher than 1.500.
- OD Control (-) is lower than 0.300.

### ◆ Results Interpretation:

The Cut Off value is 0.3.

Samples with an OD  $\geq$  than Cut Off Value are positives and samples with OD lower than Cut Off value are negatives.

The titre of the sample will be the last dilution of the sample showing an OD  $\geq$  than the Cut Off value.

## IX. ADAPTACION DEL PROTOCOLO PARA LA TITULACION A DILUCION UNICA / TITER APPROXIMATION USING A SINGLE DILUTION OF SERA SAMPLES:

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### English

Some times it is enough to have an approximate titter for the samples. In these cases, it is not needed to run 8 dilutions per sera and using a single dilution of 1/200 you may know the approximate titter.

This adjustment is based on the existing correlation between the OD value of a sample and the antibody concentration.

The procedure, to apply this interpretation, is the same as the indicated for screening. The only and important difference is the preparation of the samples. **The dilution must be 1/200** and each sample must be run in duplicate. The result of the assay will be the average between both OD values.

### Results interpretation:

Each sample will have a value of OD which will be the average between two OD values. With this average, a S/P relation must be calculated:


$$X = S/P = \frac{\text{Average of OD values of the sample}}{\text{Average of positive control OD values}}$$

$$X = S/P = \frac{\text{Media de DO de la muestra (1/200)}}{\text{Media de DO del Control positivo}}$$

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
The titre of positive sera will be determined applying the following formula:

$$y = 95,88e^{4,1x}$$

Where x is the s/p value of each sample at dilution 1/200 and "e" is the natural logarithm base (2.71828).

*It is important to remark that this is just an approximation of the titre. If there are some doubts about the results or an exact titre is had being required, we recommend to titre the samples using different dilutions.*

Desarrollado y fabricado en España por:  
Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.  
C/ Hnos. García Noblejas, 39  
28037- MADRID (SPAIN)  
Tlf: +34 91368.05.01/04  
Fax: +34 91 408.75.98  
E-mail: [ingenasa@ingenasa.com](mailto:ingenasa@ingenasa.com)  
[www.ingenasa.es](http://www.ingenasa.es)

Distributed in by: