AGROLABO
User Manual RapidVet-H Feline

AGROLABO SpA Diagnostic Division
Via Masero 59, 10010 Scarmagno (TO) - Italy
Tel. +39 0125 731111 - Fax +39 0125 731190
www.agrolabo.it - agrolabo@agrolabo.it

RapidVet-H Feline

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Feline Blood Group Determination

DESCRIPTION AND INTENDED USE

As the practice of veterinary transfusion medicine has grown significantly in recent years, the importance of identifying blood groups in cats has increased. In particular, the demand for identifying blood groups is on the rise, because only by predetermining the blood type of a blood transfusion recipient can potentially fatal transfusion mistakes be avoided. One blood group system consisting of two antigens expressed either alone or in combination has been described in cats: Type A, Type B and Type AB. The antigens are unrelated to human A B O antigen and are defined by feline alloimmune sera. Blood group incidence varies among breeds. Most cats possess the A antigen, and about one-third of those have naturally occurring, low-titered, anti-B antibody. Type B cats all have a naturally occurring, highly titered anti-A antibody. Recent studies have shown that the percentage of cats with the B antigen varies from 0.3% to 59% depending on the breed. Type AB cats are rare and since such cats have both A and B antigens on the erythrocyte membrane, they do not have or develop anti-A or anti-B antibodies.

Blood typing of cats is important in veterinary medical practice to prevent transfusion reactions in cats with A or B erythrocytes. Cats with B erythrocytes exhibit an immediate and catastrophic systemic anaphylactic reaction (hypotension, bradycardia, apnea, urination, defecation, vomiting, and severe neurological depression) and hemolytic signs (hemoglobinemia and hemoglobinuria) when transfused with Type A blood because of their natural high-titered anti-A antibody. Those cats with A erythrocytes and natural low-titered anti-B antibody will exhibit only a mild reaction when transfused with the B blood, but even this can make a difference in recovery rates in a medical situation since the transfused erythrocytes have a short life span. Other cats with A erythrocytes will not exhibit a reaction when first transfused with Type B blood but will, as a result, develop moderate titers of anti-B antibody that will result in a serious reaction upon a subsequent incompatible transfusion. Type AB cats will not exhibit transfusion reaction based on type A or B incompatibility and can receive both type A and B blood.

Blood group determination in cats is also important in making breeding decisions and in understanding medical problems in kittens. Neonatal isoerythrolysis can occur when there is blood group incompatibility between maternal and fetal blood. Because of the naturally occurring highly titered anti-A antibodies in Type B cats, neonatal isoerythrolysis can occur in Type A kittens resulting from a mating of a Type B queen with a Type A male. The maternal anti-type A antibodies occur in the colostrum where they can be absorbed by the newborn kitten, and consequently, destroy their erythrocytes. Clinically, the kittens can seem normal at birth, but develop signs after nursing, fade and die within a few days. Determining the blood groups of the queen and the male prior to mating, coupled with appropriate genetic counseling, can minimize neonatal isoerythrolysis.

RapidVet-H Feline is intended for use to classify cats as blood group Type A, Type B, or Type AB. The antisera lyophilized on a Test Card are reconstituted and well mixed with whole blood from the patient. All Type A erythrocytes react with their specific antiserum causing agglutination; all Type B erythrocytes react similarly. All Type AB erythrocytes react with both antisera causing agglutination in all cases. The results are visually identified. The characteristics of the agglutination in the "Type A" well and in the "Type B" well differs significantly because of the different nature of the antisera used.

TEST PRINCIPLE

The RapidVet-H Feline assay is based on the agglutination reaction that occurs when an erythrocyte which contains either a Type A, Type B or a Type AB antigen on its surface membrane interacts with a lyophilized antiserum specific to the particular antigen. Type A erythrocytes are characterized by the NeuGc₂GD₃ form glycolipid antigen on its surface membrane. RapidVet-H Feline uses a murine monoclonal antibody proven specific to this antigen lyophilized on the test card. The antibody molecule gives it the ability to bind and agglutinate antigens specific to Type A blood. Type B erythrocytes are characterized by the NeuAc₂GD₃ form of neuraminic acid present in the ganglioside and by the lack of the NeuGc which is present on Type A erythrocytes. The binding specificity of this form with the lectin from *Triticum vulgare* has been established. RapidVet-H Feline uses the lectin from *Triticum vulgare* to detect the presence of Type B blood.

WARNING: A certain number of feline patients exhibit auto-agglutination of varying degrees due to serum factors. Rapidvet-H Feline provides a reagent-free well to be used for the recognition of these patients ("Auto-Agglutination Saline Screen"). If a patient exhibits auto-agglutination, it will not be possible to definitively type this patient without separating the serum and serially washing the remaining red cells before performing the test. Red blood cell washing is carried out by centrifuging the sample at 2000 rpm for 5 minutes. After centrifugation, the supernatant is collected and the same volume of physiologic solution is added. Then, the operator must suspend the sample again and check if agglutination still appears. If so, he has to repeat the washing procedure again.

KIT COMPONENTS

- Agglutination Test Cards. Each card has 3 visually defined wells identified as "Auto-Agglutination Saline Screen", "Patient test Type A" and "Patient Test Type B". The cards are packaged individually in sealed polyethylene sleeves each containing a desiccant bag. Store at room temperature (+20-25 °C).
- Diluent bottle. The clear plastic bottle contains 0.02 mol/L phosphate buffered saline (PBS) at pH 7.4. The dropping tip dispenses 50 µL. Store at room temperature (+20-25 °C).
- 1 bag with pipettes and stirrers in a polyethylene bag: 1 pipette and 3 stirrers for each test.

Materials required but not provided: none.

Reagent preparation: none.

STORAGE AND SHELF LIFE

1. The agglutination Test Cards are stable at room temperature (+20-25 °C) for a period of 18 months from date of manufacture. The expiry date is printed on the back of each card.
2. The diluent is stable for 24 months from the date of manufacture. The validity of the reagent is indicated on the label.

Note: Each RapidVet-H Feline kit has an expiry date which is that of the component with the shortest validity. Other components may have a longer expiry date but their use with components of other kits should not be performed.

SAMPLE

Use feline whole blood with anticoagulant (EDTA).

Samples can be stored in the refrigerator (+2-8 °C). It is recommended to perform the test within max 5 days after drawing.

TEST PROCEDURE

1. Draw minimum 0.4 ml blood from the patient into a syringe containing EDTA as an anticoagulant. The assay requires only 150 µl of whole blood.
2. Remove the Test Card from its plastic sleeve, write the name of the cat and test date on the card. Place the Test Card on a flat surface.
3. Add 1 drop of diluent (50 µl) from the dropping bottle into the well marked as "Auto-Agglutination Saline Screen".
4. Aspirate a small amount of patient sample into the pipette and release 1 drop (50 µl) into the well marked as "Auto-Agglutination Saline Screen". This well does not contain reagents. Using a stirrer, distribute and mix the materials inside the well for about 10 seconds. (See "Note" for correct use of the pipette). A small percentage of ill cats and of healthy cats auto-agglutinate. If no agglutination is observed in the "Auto-Agglutination Saline Screen well", continue the procedure. If auto-agglutination is observed, wash the cells as indicated above (See "Warning") before proceeding.
5. Add 1 drop of diluent (50 µl) from the dropping bottle into each of the remaining wells. The diluent is required for reconstitution of the lyophilized reagents.
6. Aspirate a small amount of patient sample into a pipette and release 1 drop (50 µl) into the well marked as "Patient Test Type A", and 1 drop into the well marked as "Patient Test Type B". Using a new stirrer for each well, distribute and mix the material by gentle pressure for about 10 seconds.
7. Add a second drop of diluent in the well marked as "Patient test Type A". Do not stir with a stirrer.
8. Gently turn the card to allow the reagents on the wells to mix, taking care to avoid cross contamination between the wells "Patient test Type A" and "Patient test Type B".
9. Set the card at a 20-30° angle to allow excess blood to run to the bottom of the wells and read the results. Check for agglutination in the "Patient Test" wells.

NOTE: Use of the pipette: Hold the plastic tube between your fingers with the closed and widest side and apply pressure. Hold the specimen tube vertically and place the open end of the plastic tube on the surface of the specimen. Reduce finger pressure to draw up the sample. Hold the pipette in a vertical position directly over the well to which the sample is to be delivered. Squeeze gently and allow one free drop to fall into the well (50 µl). The pipette is designed to expel slightly in excess of 50 µl to compensate for a small amount of specimen retained by the stirrer. Use each pipette only once, then discard. Under no circumstances should the pipette be used more than once as cross-contamination can occur causing inaccurate test results.

INTERPRETATION OF THE RESULTS

If the assay has been performed correctly, visible agglutination will be present in at least one of the wells marked "Patient Test".

- If the patient sample shows agglutination in the well marked as "Type A", the cat tested is blood group A.
- If the patient sample shows agglutination in the well marked as "Type B", the cat tested is blood group B.
- If the patient sample shows agglutination in both patient wells, the cat tested is blood group AB.

Any fine, granular appearance developing after 2 minutes should be disregarded in determining the results. The kind of agglutination in the Type B well is different from that in the Type A well. Agglutination in the Type B well usually includes a small number of large, amoebic globs. Agglutination in the Type A well usually is in the form of a large number of discrete, small aggregations, each like the head of a small pin.

If the patient is very anemic, and if the patient is Type A, there may be insufficient red blood cells to obtain a significant agglutination. If the patient has a low erythrocyte volume, it might be helpful to run test without adding diluent in the patient wells (Type A and Type B).

LIMITATIONS OF THE PROCEDURE

1. To obtain reliable results it is essential that correct procedure is followed.
2. Always use a new dispensing pipette for each specimen and a new stirrer for each well. Reusing any device will cause cross-contamination and unreliable results.
3. The stability of the individual components of the kit varies. Store the components as indicated on the labels. Do not use any component after the expiration date. Use of expired materials may cause unreliable results.
4. The diluent is provided in a vial with dropper to minimize any bacterial or other contamination. Diluent from other sources in the laboratory should not be used.
5. Always draw (minimum 0.4 ml blood) with a syringe or tube containing EDTA. Less blood will cause a too high concentration of EDTA in the specimen to be tested.

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RapidVet is a trademark of dmslaboratories inc.
Agrolabo SpA