
Professional Veterinary Diagnostic Test

For Veterinary use only
For in vitro use only
Manufactured by Agrolabo S.p.A.


TEST PROCEDURE

1. Remove the test device from the package and place it on a flat surface. Keeping the **capillary pipette** in a vertical position, add 1 drop of whole blood with anticoagulant (20 µl) into the tube containing the diluent. Stir gently with the capillary pipette

2. Keeping a new **disposable pipette** in a vertical position add 9 drops of the diluted sample to the central window

3. Read the results after 5 minutes (not later than 10 minutes)

INTERPRETATION OF THE RESULTS

- If the blood group is Type A, a vertical band in “Type A” window and a horizontal band in “Control” window will appear.
- If the blood group is Type B, a vertical band in “Type B” window and a horizontal band in “Control” window will appear.
- If the blood group is Type AB, a vertical band in “Type A” window, a vertical band in “Type B” window and a horizontal band in “Control” window will appear.

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
Rev. R5 - 01/06/22

TEST PRINCIPLE

The immunochromatographic test has a three-sided layout: Side A (Type A), Side B (Type B) and Side C (Control). Specific reagents for the detection of feline erythrocytes are placed on the three membranes. The sample, diluted feline whole blood, added to the central window of the device migrates along the three membranes. The addition of the blood causes, depending on the animal's blood group and therefore on the exposed antigens, the specific reagents for group A or B bind to the erythrocytes entirely covering their surface. If the blood is Type A the erythrocytes migrate along the membrane and bind where an antibody for group A was placed. If the blood is Type B, the erythrocytes migrate along the membrane and bind with a antibody for group B. In any case, for both blood group A or B, the sample binds to the feline anti-erythrocyte antibody adhering to the membrane at the level of the window C (control), creating a red band showing the correct test performance.

KIT COMPONENTS

- Disposable capillary pipettes
- Disposable pipettes
- Tubes with diluent
- RapidVet-H IC devices

STORAGE AND SHELF LIFE

Store at room temperature (+18-25°C). If correctly stored, the product remains stable until expiry date. Avoid the kit being exposed to excessive heat or very low temperatures.

SAMPLES

Whole blood with anticoagulant (EDTA, heparin), fresh or stored in the fridge for no longer than 5 days, must be used. Haemolytic samples can interfere with the test results.

PRECAUTIONS

Test procedure must be followed carefully. Open device packaging immediately prior to use. For veterinary use only. Do not use test components after expiry date. The test device should be placed horizontally on a flat surface while testing. Use a new disposable pipette for each sample. Always take min 0.4 ml of blood to be added into the test tube with anticoagulant. Smaller blood quantities could lead to a high concentration of anticoagulant in the sample to be tested. Use samples avoiding clots and haemolysis, in order to obtain correct results. Use the diluted sample immediately after dilution.

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RapidVet-H IC Feline - Immunochromatographic test for blood group determination in cat (A, B, AB)