Professional Veterinary Diagnostic Test
BioPronix product line
For Veterinary use only
For in vitro use only
Manufactured by Agrolabo S.p.A.

TEST PROCEDURE

1. Bring all test samples and kit components to room temperature at least 10 minutes before the test is to be performed.
2. Remove the device from the pouch and place it on a clean, flat surface.
3. Keeping the **capillary pipette** in a vertical position, add 1 drop of whole blood with anticoagulant inside the tube containing the diluent. Stir gently with the capillary pipette.
4. With a new **disposable pipette**, keeping it vertically, add 4 drops of the diluted sample into the S window of the device.
5. Read the results after 5 minutes (no more than 10 minutes).

QUICK GUIDE


1. Remove the device from the wrapper and place it on a clean, flat surface. Keeping the **capillary pipette** in a vertical position, add 1 drop of whole blood with anticoagulant inside the tube containing the diluent. Stir gently with the capillary pipette.
2. With a new **disposable pipette**, keeping it vertically, add 4 drops of the diluted sample into the S window of the device.
3. Read the results after 5 minutes (no more than 10 minutes)
4. Interpretation of results

INTERPRETATION OF THE RESULTS

- If the sample blood is type A, 2 lines will appear within the central window: one line A and one line C (control).
- If the sample blood is type B, 2 lines will appear within the central window: one line B and one line C (control).
- If the sample blood is type AB, 3 lines will appear within the central window: one line A, one line B, and one line C (control).
- Invalid Result: the test is considered invalid if no line appears in the central window at line C (control), even if line A, B or both are present.

WARNING Open the test package immediately before use. Place the device horizontally for testing. Keep the pipette in a vertical position when adding the sample to the device.

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
Professional Veterinary Diagnostic Test
BioPronix product line
 For Veterinary use only
 For in vitro use only
 Manufactured by Agrolabo S.p.A.

Feline AB IC - Rapid single-membrane immunochromatographic test for determination of cat blood type (A, B, AB) with new monoclonal antibodies

USER MANUAL

Rev. R0 -24/04/24

TEST PRINCIPLE

Feline AB IC is an immunochromatographic test for the rapid determination of blood type in cat whole blood.

Cat erythrocyte-specific reagents, capable of detecting the presence of exposed A and B antigens on red blood cells, are immobilized on a single strip.

The sample, diluted cat whole blood, added to the S window of the device, migrates along the membrane and binds to the monoclonal antibodies specific for A or B, depending on the animal's blood type.

All results are displayed within the central window.

If the blood is type A, the erythrocytes migrate across the membrane and bind where the antibody for group A is present (line A). If the blood is type B, the erythrocytes migrate across the membrane and bind where the antibody for group B is present (B line). If the blood is type AB, the erythrocytes migrate across the membrane forming both A and B lines.

In any case, with both type A and type B blood, the feline sample binds to the anti-erythrocyte antibody adhered at line C (control), creating a red line that represents the test's functioning control.

KIT COMPONENTS

- Feline AB IC devices
- Tubes with diluent
- Disposable capillary pipettes
- Disposable pipettes

STORAGE AND SHELF LIFE

Store at room temperature (+18-25°C). If correctly stored, the product remains stable until expiry date. Avoid the kit being exposed to excessive heat or very low temperatures. Do not freeze.

SAMPLES

Perform the test using whole blood with anticoagulant (EDTA, heparin) fresh or stored in a refrigerator for no more than 5 days. Inaccurate results are possible in the presence of autoimmune hemolytic anemia or due to improper sample storage.

PRECAUTIONS

- Test procedure must be followed carefully
- The test device should be placed horizontally on a flat surface while testing
- Use the diluted sample immediately after dilution
- Use a new pipette for each sample
- The presence of moisture may decrease the stability of the reagents, so open the device package only at the time of use.
- Do not use the components after the expiration date and do not reuse the components in the kit or from other devices.
- Do not touch the membrane inside the central results window.
- For veterinary use only.