AGROLABO
FilarCHECK Ag ELISA

FilarCHECK Ag ELISA

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
Tabella 1

FilarCHECK Ag ELISA

Kit ELISA for the detection of Dirofilaria immitis antigen in canine serum or plasma samples

TEST PRINCIPLE

The test is performed as a direct sandwich enzyme linked immunosorbent assay (direct sandwich ELISA). A specific monoclonal antibody (Mab) against the antigen of *Dirofilaria immitis* is fixed on the wells of the plate. The samples to be analyzed are distributed in each well. If the sample contains the specific antigen of *Dirofilaria immitis*, the antigen will bind to the Mab adsorbed, forming the Ag-Ab complex. By adding a second Mab anti-*Dirofilaria immitis* (conjugated with horseradish peroxidase), this will bind to the Ag-Ab complex. All non-fixed material is removed by washing. A specific substrate is then added to the wells which causes a colorimetric reaction.

TEST COMPONENTS

- 1 x 48 or 96 wells plate (in strip of 8 breakable wells)
- 1 dropper of positive control - ready for use
- 1 dropper of negative control - ready for use
- 1 dropper of conjugate - ready for use
- 1 bottle of washing solution - 10X
- 1 dropper of substrate (TMB) - ready for use
- 1 dropper of stop solution - ready for use

MATERIALS REQUIRED BUT NOT SUPPLIED

Test tubes for serum dilutions, precision pipettes and tips, spectrophotometer at a wavelength of 450 nm.

PRECAUTIONS

1. **Read carefully and follow the instruction for use**
2. Do not use reagents and instructions from other kits
3. Bring all reagents to room temperature (+20-25 °C) 1 hour before running the test
4. Perform the test at room temperature and keep it constant throughout the execution
5. Do not use kit and reagents after the expiration date
6. Avoid thermal stress: use all reagents for the time strictly necessary for the test and store them at +2-8 °C as soon as possible
7. Take out only the strips needed to perform the test, reseal the remaining unused ones in the appropriate bag with desiccant, seal and store at +2-8 °C
8. Prepare samples according to instructions: if possible, avoid using lipemic and hemolyzed samples
9. Always use the positive and negative control in each analysis
10. Mix reagents before use
11. Respect the times of each incubation and the order of deposition in the plate
12. Use a new tip for each sample/reagent to be analyzed, avoiding contamination
13. TMB substrate is very sensitive to light and contamination. Do not insert the tip directly into the vial. Handle with care
14. The stop solution is a strong diluted acid. Handle with care
15. All components should be stored at +2-8 °C
16. Do not smoke, eat or drink where specimens or kit reagents are being handled

PREPARATION OF REAGENTS

Wash buffer solution:

Dilute one part of the wash buffer supplied into 9 parts of distilled or deionized water (ie. 100 ml of concentrated wash buffer with 900 ml of water). Once diluted, the wash buffer solution must be stored at +2-8°C.

Note: the 10X wash buffer solution may have some precipitates at the bottom of the bottle that can be solubilized by bringing the solution to room temperature or placing it at +37°C for about 15 minutes.

WASHING STEPS

The washing steps may be carried out using an automated plate washer or a multichannel pipetting device suitable for dispensing 300 µl in each well. After the incubation, the washing steps should be carried out according to the following instructions:

- throw out the content of the plate by vigorously turning the plate over to avoid the possible mixture of the content from one well to another
- dispense a volume of 300 µl of diluted wash buffer into each well
- shake the plate gently, avoiding contamination between wells
- turn over the plate and tap the wells vigorously to remove the wash buffer
- repeat the process as many times as is indicated by the kit's instructions
- prior to emptying out the plate's contents during the final washing step, ensure that the next reagent to be added to the plate is ready to use
- the wells should not be left dry for any longer than strictly necessary in order to dispense the reagent of the following step
- after the last washing step, blot the empty wells face down on clean paper towels or absorbent paper to remove residual wash buffer

PREPARATION OF SAMPLES

Canine serum or plasma samples.

Dilution of the samples is not necessary. Directly add 30 µl of sample to each well.

TEST PROCEDURE

1. Bring all reagents to room temperature 1 hour before use.
2. Dispense 30 µl of positive control and 30 µl of negative control in two different wells. Add 30 µl of each sample to the remaining wells of the plate. Dispense 70 µl of conjugate into each well, following the order in which the samples and controls were dispensed. **Delicately shake the plate, seal and incubate at room temperature for 10 minutes.**
3. Wash cycle 5 times according to the procedure previously described.
4. Add 70 µl of substrate in each well and delicately shake the plate. **Seal the plate and incubate at room temperature for 10 minutes in the dark.**
5. Proceed with visual or spectrophotometric reading as explained in points A and B:

- A. Visual reading:** at the end of the incubation with substrate, read visually. Evaluating blue colour in the wells.
- B. Spectrophotometric reading:** at the end of the incubation with substrate, add 70 µl of stop solution to each well. Read the results with a spectrophotometer at a wavelength of 450 nm.

READING AND INTERPRETATION OF THE RESULTS

A. VISUAL READING - Without adding stop solution

The positive control must be blue/light blue and clearly distinguishable from the negative control, which should be colourless.

Interpretation of the results

POSITIVE samples: a blue colour will develop. Colour intensity is proportional to the quantity of antigen present in the sample.

NEGATIVE samples: no colour will develop.


B. SPECTROPHOTOMETRIC READING - Add stop solution

The reading of the absorbance values (OD) must be carried out at a wavelength of 450 nm. If you run the samples or controls in duplicate wells, the OD values for the samples or controls will be the mean of both wells.

Interpretation of the results

Samples with OD values ≥ 0.100 must to be considered positive. (See Table 1).

Table 1

AGROLABO
FilarCHECK Ag ELISA

Manufactured by
Agrolabo SpA

Head Office, Laboratories
and Production Centre
Diagnostic Division
Via Masero 59
10100 Scarmagno (TO)
Italy
Tel +39 0125 731111
E-mail agrolabo@agrolabo.it
www.agrolabo.it