## INGEZIM IBR 2.0

Prod Ref: 12.BHV.K1

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## I. TECHNICAL BASIS

The kit has been designed to detect specific antibodies against BoHV1 (Bovine Herpesvirus 1) in affected cattle of infectious bovine rhinotracheitis (IBR). This Kit is based on an indirect enzymatic immunoassay. We make a brief description of the technique below:

The antigen is fixed in a solid support (polystyrene plate). The sample is added to each well. After incubation we add a monoclonal antibody (peroxidase conjugated) specific to bovine immunoglobulins. If the sample contains antibodies against the virus, the

conjugate will bind to the plate, whereas if it does not contain specific antibodies no binding of the conjugate will be appreciated.

After washing the plate to eliminate all non-fixed material, we can detect the presence or absence of labelled conjugate by adding the substrate which, in presence of the peroxidase will develop a colorimetric reaction.

The assay is suitable for use both as sera and milk, whether individual or tank

## II. PRECAUTIONS AND WARNINGS FOR USERS:

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature (20°-25°C) prior to use.
3. Do not mix instructions or reagents from different kits.
4. Avoid any contamination of the kits reagents.
5. Do not use components after expiration dates and do not mix components from different lots.
6. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled.
7. Do not pipette by mouth.
8. Use a new tip for each samples sample.
9. Stop solution is a strong acid solution that must be used with precaution. In case of accidental contact with skin, wash gently with water.
10. For each utilisation of the Kit, control positive and negative samples must be tested in a systematic way.
11. Substrate must be handle with care, it is very sensible to light and contamination.

## III. STORAGE OF COMPONENTS

All reagents and plates must be stored between +2°C and +8°C. **Once opened**, control are stable for one month. In case that they are not going to be used in this period, we recommend to make aliquots and store them at -20°C.

## IV. INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µl on each well.

After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turn over of the plate to avoid the possible mixture of the content from one well to another.
- Dispense a volume of 300 µl of washing solution on each well.

Shake delicately the plate, avoiding the contamination between wells.

- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit.
- Prior to empty the content of the last washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry more time than strictly needed.
- After the last step of washing shake the plate turned over an absorbent filter paper.

## V. PREPARATION OF SAMPLES:

### Sera:

Samples must be diluted 1/5 in sample diluent DE01-01. The dilution can be done in the same well by first adding 80 µl of diluent and after that 20µl of serum.

### Milk (individual and tank)

Milk must be tested undiluted. Fresh, refrigerated or previously frozen milk may be tested. To eliminate interference from lipids, partially skimmed milk should be either by centrifugation (15 min at 2000xg) or leaving overnight at 4°C to observe the formation of a lipid layer on the surface. Sample should be collected under the layer of lipids and that will be the one used in the assay. To obtain whey, freeze/thaw the milk sample.

## VI. PREPARATION OF REAGENTS

- **Washing solution:** Dilute one part of the concentrate washing solution provided in the kit with 24 parts of distilled or deionized water (40 ml of concentrate solution and 960 ml of water). Once this solution ready, it remains stable between +2°C and +8°C.

### Controls

Controls will be diluted 1/5 in diluent DE01-01 and will be used at this dilution whether the test is done with sera or with milk.

Once opened we recommend freezing the control at -20°C if they are not going to be used in a short period of time.

### Preparation of the conjugate: to make immediately before use:

Dilute the needed quantity of conjugate provided in the Kit 1/100 with conjugate diluent (DE08-01)

**For serum samples (conjugate with yellow cup)**  
**For milk samples (conjugate with green cup)**

**ATTENTION:** Shake the solution before use. Prepare only the quantity needed for each time because the remainder volume has to be rejected.

## VII. TEST PROCEDURE

1. All reagents (except concentrated conjugate) must be equilibrated at room temperature before use
2. **Sample addition (serum or milk):**
3. Dispense 100 µl/well of each sample prepared as indicated in section V. Add last dilute controls as indicated in section VI. It's recommended to evaluate samples and controls in duplicate. **Incubate 45 min at 37°C.**
4. Wash 4 times following the described procedure.
5. Add 100 µl of conjugate (prepared as described in section VI) to each well. **Incubate the plate for 45 at 37°C.**
6. Wash 5 times following the described procedure.
7. Add 100 µl of substrate solution, to each well. Incubate: **SERUM: 15 min and MILK: 5 min**
8. Incubations must be performed at room temperature in the dark. Count the time from the addition of the solution to the first well. It's recommended using a multichannel pipette to expedite this process as much as possible
9. Add 100 µl of stop solution to each well. **CAUTION:** The stop solution must be dispensed in the same order as the substrate solution was added.
10. Read the OD with a spectrophotometer at 450 nm within 5 min after the addition of stop solution.

## VIII. READING AND RESULT INTERPRETATION

The reading must be done with a spectrophotometer at **450 nm**. In case that samples had been tested in duplicate, it must be considered the average of the two OD values obtained. Similarly, the arithmetic mean of the values obtained in the two wells for positive control and two wells for the negative control should be calculated.

### A. Validation criteria:

The test can be considered valid if:

#### SERA:

- The OD of Positive Control (1/5 dilution) > 0.8
- The OD of Negative Control (1/5 dilution) < 0.3

#### MILK:

- The OD of Positive Control (1/5 dilution) > 1.0
- The OD of Negative Control (1/5 dilution) < 0.40

### B. Cut off calculation:

C.

#### SERUM:

First, calculate the S/P ratio (sample/positive) using the following formula

$$S/P = (S-NC) / (PC-NC)$$

Where S is the OD sample value and NC and PC are the OD values of negative and positive controls respectively.

- **Positive/Negative Cut Off**

**The cut-off value is 0.18.**

Samples showing S/P values equal or higher than the cut off are considered positive and samples showing S/P values lower than the cut off are considered negative.

#### MILK:

- **Positive/Negative Cut Off**

**Mean of OD450 of Negative control + 0.1.**

Samples showing OD values equal or higher than the cut off are considered positive and samples showing OD values lower than the cut off are considered negative.

Developed and manufactured in Spain by:

INMUNOLOGIA Y GENETICA APLICADA, S.A.  
C/ Hnos. García Noblejas, 39  
28037- MADRID (SPAIN)  
Tif: +34 91368.05.01/04  
Fax: +34 91 408.75.98  
E-mail: [ingenasa@ingenasa.com](mailto:ingenasa@ingenasa.com)  
[www.ingenasa.es](http://www.ingenasa.es)

IT-73840  
IT-73780

ISO 14001:2015  
9191.INGE

ISO 9001:2015  
9175.ING2

Distributed in

by: