
biopronix
by **agrolabo** S.p.A.

DermaKit

**In-house screening test for the rapid identification
of dermatophyte infections in dogs, cats and horses.**

AGROLABO SpA Diagnostic Division

Via Masero 59, 10010 Scarmagno (TO) - Italy
Tel. +39 0125 731111 - Fax +39 0125 731190
www.agrolabo.it - agrolabo@agrolabo.it

DermaKit

In-house screening test for the rapid identification of dermatophyte infections

INTRODUCTION

Dermatophytosis, often called “ringworm”, is a cutaneous infection caused by different genera of fungi collectively called “dermatophytes”, meaning “plants that live on the skin”. The main fungi responsible for dermatophytosis in domestic animals are *Microsporum*, *Trichophyton* and *Epidermophyton*.

These pathogenic fungi are found worldwide and are the most frequent cause of cutaneous mycosis in companion animals, as well as in humans. Dermatophytosis is important as it is a highly contagious zoonotic skin disease. Children, the elderly and the immunocompromised are special at-risk populations, but anyone in frequent contact with infected pets risks contracting the disease.

INTENDED USE

DermaKit is intended for use as part of a differential diagnosis of dermatological conditions in dogs, cats and horses.

DESCRIPTION

For the isolation and general identification of dermatophytes, Taplin and colleagues formulated a culture medium known as a Dermatophyte Test Medium (DTM).

DermaKit is a Dermatophyte Selective Medium (DSM) based on the original Dermatophyte Test Medium. The newly formulated DSM contains:

- specific nutrients that promote the growth of the dermatophytes.
- selective antibiotics that inhibit the growth of the non-pathogenic saprophytic mycetes and of bacteria.
- a pH indicator that implies dermatophyte presence through a colour change to red as a result of the *Microsporum*, *Trichophyton* and *Epidermophyton* genera producing alkaline metabolites.

KIT CONTENTS

- 12 tubes, Labels and Work station

STORAGE

If the test is stored at +2-8°C, stability is guaranteed for 36 months (the expiry date is indicated on the external label of the kit). If the test is stored at room temperature, stability is guaranteed for about 20 months.

TEST PROCEDURE

- Inspect the site of sampling to ensure that it has not been treated with medications that could affect results.
- Bacterial or saprophytic contamination can be reduced by gently cleaning the area from which the sample is taken with 70% alcohol.
- Using a sterile scalpel blade, take scales and/or hair from the border of lesion.
- Cats are frequently asymptomatic carriers of dermatophytes that can infect other cats or man. In order to test these animals for the presence of dermatophytes, use a new (or sterilized) toothbrush to brush all the body or possible suspected areas and place the hair thus obtained on the culture medium.
- Place the sample on the surface of the culture medium without cutting into it.
- Close the tube without completely tightening the lid, to avoid the formation of moisture that could alter the result.
- Incubate at room temperature (+22-25°C) for 14 days. Check daily for a change in the medium's colour and for colony growth

INTERPRETATION OF THE RESULTS

The medium contains a phenol red indicator, which changes from yellow to red when dermatophyte fungi are present in the patient sample. In order to survive and grow, dermatophyte pathogens immediately use the protein source in the medium, determining a rapid alkalization of the medium and consequently provoking the colour change. In such cases, the medium will change colour before colony growth is visible. Contaminants, on the other hand, generally use the carbohydrate source first and only after exhausting carbohydrates they use proteins. Thus, a colour change is not produced until after all of the carbohydrate source has been used. Contaminant colonies, in fact, grow and are visible on the media before the test changes colour.

Periodically examine the reaction substrate in each tube. Any change in colour from yellow to red, before colony growth, even in only a small region of the reaction substrate, is interpreted as positive for the presence of dermatophytes. However, the degree of red will normally intensify, and the area of colour will normally spread with time.

Special nutrients have been added to the **DermaKit** medium to help promote the growth of dermatophytes. Statistically positive results will be visible in approximately 72 hours. Depending on various factors, however, including the type of microorganism, the amount of dermatophyte material sampled and the stage of the disease, a longer incubation period may be required.

MACROSCOPIC FEATURES OF THE MAIN DERMATOPHYTE COLONIES

| Genus Tricophyton | Colony description |
|--|---|
| T. mentagrophytes | grainy, floury, feathery and cottony of white colour, then yellow and pink colour. Reverse side: pink brown |
| T. rubrum | flat or floury and cottony velvety heaps of white cream colour, then dark pink. Reverse side: purple red |
| T. tonsurans | powdery and velvety, first flat then in heaps, coloured from cream to nut or from yellow to pink. Reverse side: from yellow to mahogany red |
| T. violaceum | heaped and verrucous, cream and then purple or violet colour |
| T. equinum | white, filmy with yellow pigment on border then cream and nut colour. Reverse side: first yellow then pink and red brown |
| T. schoenleinii T. ferrugineum T. concentricum | irregular, hard, white then nut colour, of waxy, powdery and feathery aspect |

| Genus Microsporium | Colony description |
|---------------------------|--|
| M. audouinii | velvety, flat, brownish colour. Reverse side: light orange colour |
| M. gypseum | powdery with fringed borders of cinnamon colour |
| M. canis | cottony, white and yellow on border. Reverse side: yellow then orange-brown colour |

| Genus Epidermophyton | Colony description |
|-----------------------------|---|
| E. floccosum | downy white then velvety and then greenish. Reverse side: dark yellow |

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Biopronix Product Line by Agrolabo S.p.A.