**biopronix**  
by **agrolabo** s.p.a.

## **MalasseziaKit**

**In-house screening test for the detection  
of Malassezia infection in cats and dogs**

# MalasseziaKit

## *Test for the detection of Malassezia*

### **INTRODUCTION**

Malassezia is a yeast which is commonly found in the mucous membranes and the skin of dogs and cats.

This micro-organism is a part of the normal microflora. However it can become pathogenic under certain conditions and cause otitis and severe itchy dermatitis.

Predisposing factors can cause a pathological proliferation of Malassezia. Malassezia dermatitis and otitis often worsen atopic dermatitis and can be associated with infection by Staphylococcus. Moreover, endocrine and other metabolic disorders can cause seborrhea, which can lead to an abnormal overgrowth of Malassezia.

Pathogenicity of Malassezia is very complex due to the fact that it is often accompanied or caused by the diseases the differential diagnosis is based on. Therefore, the diagnosis should be confirmed in order to start targeted treatments.

### **DIAGNOSTIC METHOD**

Diagnosis of Malassezia dermatitis or otitis is made by clinical examination, cytology and fungal culture.

The cytological examination is carried out on samples taken from the ear canal or the skin surface. These are applied to a slide and examined microscopically, with special stains.

The minimum number of yeasts indicating the possibility of a real infection is not standardized. This technique has a very low sensitivity which leads to false negative results.

The fungal culture is used to confirm the results of cytology or the hypothesis of infection, if the cytological examination is negative. Even in this case the minimum number of colonies indicating the real infection is not standardized.

## DESCRIPTION

*Malassezia pachydermatis* is the most widespread species in dogs and cats and is the only species which is generally lipophilic. Other species of *Malassezia* (*furfur*, *globosa*, *sympodialis*) have been isolated in dogs and cats, but these are lipodependent, which means they cannot be identified by normal culture media.

A specific chromogenic culture medium has been studied. This allows for the growth and the identification of all *Malassezia* species in a few days.

## KIT CONTENTS

- 12 tubes,
- Labels

## STORAGE AND SHELF LIFE

Store the test at + 2-8°C.

Stability is guaranteed until the expiry date. Store in the dark.

Bring the vial to room temperature 1 hour before use.

## TEST PROCEDURE

**Otitis:** samples of earwax from infected ears can be taken with a cotton swab. The ideal practice would be to collect ear exudates in the horizontal portion of the canal. In this way, the material collected is more abundant and the results are more representative.

**Dermatitis:** take skin scrapings from the affected areas where erythema, lesions or pruritus are present with a swab. If possible, dip the swab into saline solution.

- Inoculate the sample slithering it in the culture medium, as soon as possible after collection. It is recommended not to refrigerate the samples after the collection in order not to affect the yeast vitality.
- Place the cap on the vial, without tightening it.
- Incubate at +32°C for 5 days. Check daily for a change in the medium colour and for colony growth.

If *Malassezia pachydermatis* or other *Malassezia* species are pathogenic, the growth of pink, dark pink or violet smooth rounded colonies occurs.

## RESULT INTERPRETATION

If *Malassezia pachydermatis* is present, in most of the cases pink, dark pink or violet smooth rounded colonies grow within 48-72 hours. If more than 10 colonies are visible and clinical signs occur, this indicates that an infection is in progress.

Other *Malassezia* species can grow and be identified according to the colony colour. In most of the cases these can be found in association with *Malassezia pachydermatis*.

Neither other bacterial nor dermatophyte contamination can develop, but the growth of *Candida* spp whose colonies ranges from green to light blue is possible.

Yeast identification can be made by direct observation with the microscope. *M. pachydermatis* is typically peanut-shaped.

## INTERPRETATION OF THE RESULTS

*Malassezia pachydermatis* (smooth, from pink to violet)

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**Biopronix Product Line by Agrolabo S.p.A.**