Professional Veterinary Diagnostic Test  
**BioPronix** product line  
 For Veterinary use only  
 For in vitro use only  
 Manufactured by Agrolabo S.p.A.

## INTERPRETATION OF THE RESULTS

**Negative Results:** The test is negative if a red line appears only in window C (control line).

**CPV Positive Results:** The test is positive for CPV if two red lines appear in column 1. One line will appear in window T (test line) and one in window C (control line).

**Giardia Positive Results:** The test is positive for Giardia if two red lines appear in column 2. One line will appear in window T (test line) and one in window C (control line).

**CCoV Positive Results:** The test is positive for CCoV if two red lines appear in column 3. One line will appear in window T (test line) and one in window C (control line).

## QUICK GUIDE


1. Take a faeces sample and insert it into the tube. Shake well and place the tube in a vertical position.
2. Once the residues have settled, collect the supernatant with one of the pipettes provided in the kit and add 3-4 drops into all 3 S windows.
3. Wait 5-10 minutes to read the results, not later than 15 minutes.
4. Interpretation of the results.

### Invalid Result:

The test is invalid if no coloured band appears in window C, even if a coloured band appears in window T.

### WARNING

Open the test package immediately before use. Place the device horizontally for testing. Keep the pipette in a vertical position when adding the sample to the device.

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
Negative

Invalid

## PARVO/GIARDIA/CORONA IC - Rapid immunochromatographic test for the detection of canine Parvovirus, Giardia lamblia and canine Coronavirus in faeces samples

## USER MANUAL

Rev. R2 – 02/02/21

### TEST PRINCIPLE

Parvo/Giardia/Corona IC is a qualitative test that employs the immunochromatographic sandwich technique using two monoclonal antibodies that recognize two different epitopes of the antigens of canine Parvovirus (CPV), Giardia lamblia and canine Coronavirus (CCoV). A monoclonal antibody is conjugated with colloidal gold and a further monoclonal antibody is immobilized on the test membrane.

If the sample contains CPV, Giardia and CCoV antigens (positive sample), these antigens will bind to the colloidal gold conjugated antibody forming the Ab-Ag complex. The complex migrates along the membrane and will be captured by the second monoclonal antibody immobilized on the membrane at T window level, forming a red line (test line). The absence of the red line in the T window will indicate a negative result. The sample continues to migrate along the membrane up to window C, where a reagent able to bind the antibody conjugated with colloidal gold is immobilized, resulting in the formation of a red line (control line). The presence of the red line in window C confirms that the test was performed properly and should always appear, with both positive and negative samples.

### KIT COMPONENTS

- Disposable pipettes
- Tubes with diluent
- Swabs
- PARVO/GIARDIA/CORONA IC devices

### STORAGE AND SHELF LIFE

Store at room temperature (+18-25°C). When properly stored, the product remains stable for its full shelf life. Keep the kit away from excessive heat/cold or direct sunlight.

### PRECAUTIONS

**IMPORTANT:** open individual device packaging immediately prior to use. For veterinary use only. Do not use kit components after expiry date. The test device should be placed horizontally on a flat surface while testing. Use a new disposable pipette for each sample.

### TEST PROCEDURE

1. Take a faeces sample dipping the swab into the faeces of the dog (not more than 30 mg). Insert the sample into the appropriate tube containing diluent. Carefully shake the tube to mix the substances and obtain the correct sample dilution. Place the tube in a vertical position, in order to allow residues to precipitate. The faeces samples can be stored at +2-8°C if used within 24 hours or at -20°C for analysis to be done later.
2. Once the residues have precipitated, use a new pipette to extract the supernatant from the tube. Remove the test device from its individual wrapping, placing it on a flat surface.
3. Holding the pipette approximately 1 cm from the test device, add 3-4 drops of supernatant to all 3 windows S of the device.
4. Read the results after 5 -10 minutes and no later than 15 minutes