



ELISA KITS

ELISA kits for the detection
of the main pathologies in dogs and cats



ELISA Test

Enzyme-Linked Immunosorbent Assay Technology

Agrolabo manufactures a series of direct ELISA kits, employing monoclonal antibodies to detect the presence of a particular antigen in biological samples and indirect ELISA kits, determining the presence of specific antibodies in serum samples.

- Agrolabo's ELISA tests are highly sensitive and specific, designed to screen large numbers of samples at a time
- Reliable and accurate results
- Our ELISA tests are practical, easy to use
- Each of our ELISA tests is subjected to a careful validation examination before being offered on the market
- Additional equipment such as a plate reader, may be required when you prefer the spectrophotometric reading at 450 nm
- Our Technical Service Team may assist you in the choice of most suitable product
- Each kit is complete with all reagents ready to use
- Format 96 wells (breakable wells)
- Visual or spectrophotometric reading of the results
- 18-months shelf life from the date of production
- Storage at +2-8°C



Canine ELISA Test

One Step Test

EhrlichiaCHECK Ab ELISA

Indirect ELISA test for the detection of IgG antibodies against Ehrlichia canis
Sample serum or plasma 96 Wells
Sensitivity 93.9% - Specificity 100%

FilarCHECK Ag ELISA

Direct Sandwich ELISA test for the detection of Dirofilaria immitis antigen
Sample serum or plasma 96 Wells
Sensitivity 97.6% - Specificity 100%

LeishmaCHECK Ab ELISA

Indirect ELISA test for the detection of IgG antibodies against Leishmania infantum
Sample serum or plasma 96 Wells
Sensitivity 95.5% - Specificity 99.5%

DistemperCHECK IgG Ab ELISA

Indirect ELISA test for the detection of IgG antibodies against Canine Distemper Virus (CDV)
Sample serum or plasma 96 Wells
Sensitivity 99% - Specificity 100%

ParvoCHECK IgG Ab ELISA

Indirect ELISA test for the detection of IgG antibodies against Canine Parvovirus (CPV)
Sample serum or plasma 96 Wells
Sensitivity 95% - Specificity 98.5%

ParvoCHECK Ag ELISA

Direct Sandwich ELISA test for the detection of Canine Parvovirus (CPV) antigen
Sample faeces 96 Wells
Sensitivity 99% - Specificity 100%

Feline ELISA Test

One Step Test

FeLVCHECK Ag ELISA

Direct Sandwich ELISA test for the detection of Feline Leukemia Virus (FeLV) p27 antigen
Sample serum or plasma 96 Wells
Sensitivity 97.6% - Specificity 97.1%

FIVCHECK Ab ELISA

Indirect ELISA test for the detection of IgG antibodies against Feline Immunodeficiency Virus (FIV)
Sample serum or plasma 96 Wells
Sensitivity 100% - Specificity 100%

FCoVCHECK Ab ELISA

Indirect ELISA test for the detection of IgG antibodies against Feline Coronavirus (FCoV)
Sample serum or plasma 96 Wells
Sensitivity 93.5% - Specificity 100%

FeLVCHECK gp70 Ab ELISA

Indirect ELISA test for the detection of antibodies against gp70 of Feline Leukemia Virus (FeLV)
Sample serum or plasma 96 Wells
Sensitivity 96% - Specificity 100%





Head office, laboratories and production centre:
Via Masero 59 - 10010 Scarmagno (TO) Italy - Tel. +39 0125 731111 - agrolabo@agrolabo.it
www.agrolabo.it

COMPANY WITH CERTIFIED QUALITY
MANAGEMENT SYSTEM
= ISO 9001 =